

## CMV

## Cytomegalovirus DNA by PCR - Quantitative

<b>GA Test Code</b>	<b>3702</b>
<b>Method</b>	Quantitative Real-Time Polymerase Chain Reaction (qPCR)
<b>Specimens</b>	<p><b>Urine:</b> 10.0 (min 5.0) mL, refrigerated (7 days).</p> <p><b>CSF:</b> 1.0 (min 0.25) mL, refrigerated (7 days) or frozen.</p> <p><b>Swab (e.g. newborn saliva, or from any other site):</b> Collect sample and place entire swab in 2.0 mL saline or viral transport media in a sterile screw top tube. Do not use calcium alginate or wood shafted swab. Ship ambient up to 14 days.</p> <p><b>Whole Blood (ACD or EDTA):</b> 5.0 (min 3.0) mL, ambient (4 days), refrigerated (7 days).</p> <p><b>Plasma (ACD, EDTA, or PPT):</b> 3.0 (min 1.0) mL, separated/centrifuged within 6 hours, refrigerated or frozen (<i>do not freeze in PPT</i>). If storing longer than 24 hours, store frozen.</p> <p><b>Fluid (e.g. amniotic, peritoneal, pleural):</b> 2.0 (min 1.0) mL, ambient (4 days).</p> <p><b>Bronchial Washing:</b> 3.0 (min 1.0) mL, refrigerated (7 days).</p> <p><b>Sputum:</b> 10.0 (min 5.0) mL, refrigerated (7 days).</p> <p><b>Stool:</b> 4-8 g of feces, screw-cap container, refrigerated (7 days). Do <b>not</b> dilute the specimen or use preservatives.</p> <p><b>Other Samples:</b> Please contact GA for questions about other specimens.</p>
<b>Causes for Rejection</b>	Quantity not sufficient (QNS) for analysis; time and/or temperature instructions not followed; blood in heparin; plasma frozen in PPT; calcium alginate or wood shafted swab; no swab in tube and/or received ambient after 14 days.
<b>Reference Range</b>	Not Detected (< 500 IU/mL)
<b>Quantitative Range</b>	500 to 2.5 x 10 <sup>10</sup> CMV DNA IU/mL
<b>Turnaround Time</b>	Same or Next Day
<b>CPT Code</b>	87497

### Description

Cytomegalovirus (CMV) DNA is detected by a real-time PCR assay utilizing PCR primers directed against viral sequences found in the US17 region of the CMV genome. A patient value of less than 500 CMV DNA IU/mL indicates that the patient's viral load is below the quantitative limit of this assay, but does not indicate that the patient is not infected with CMV.

### Clinical Utility

CMV is a commonly found virus that threatens immunocompromised patients including neonates, transplant recipients, oncology patients and patients with AIDS. Commonly seen manifestations of a CMV infection include: encephalitis, retinitis, colitis, hepatitis, adrenalitis, polyradiculopathy, and esophagitis. CMV is the major viral pathogen that causes death after renal transplantation. The use of PCR has been found to detect CMV infection at a much higher rate in renal allograft cases, thus resulting in improved patient management.

Every year, 1 in 150 children is born with congenital CMV infection, resulting in possible hearing loss. Studies have shown that using a real-time PCR assay to screen newborn saliva for CMV yielded at least 97.4% sensitivity and 99.9% specificity when compared to culture. CMV-infected babies can be monitored closely for hearing loss, with support services made available as necessary.

Boppana S, et al. Saliva polymerase-chain-reaction assay for cytomegalovirus screening in newborns. *N Engl J Med* 2011; 364: 2111-2118.

Liapis, et al. CMV infection of the renal allograft is much more common than pathology indicates: a retrospective analysis of qualitative and quantitative buffy coat CMV-PCR, renal biopsy pathology and tissue CMV-PCR. *Nephrol Dial Transplant* 2003; 18:397-402.